

# Moving Posture Pro Records

All the patient information for Posture Pro is stored in a folder named PPDATA. In pre Posture Pro 8 versions this folder was usually at C:\PPV\PPDATA. In Posture Pro 8 version it's located at C:\PP8\PPDATA.

To move data from one computer to another you will copy the PPDATA folder and drop or paste it to the new location of either C:\PPV or C:\PP8. This will overwrite any existing data in the PPDATA folder of the new machine, so perform this task before adding any new records to the new version of Posture Pro.

You can also use the Backup/Restore Feature of Posture Pro

<http://www.posturepro.com/help/Backup%20or%20Restore%20Database/Backup%20or%20Restore%20Database.html>

Online URL: <https://posturepro.phpkb.cloud/article.php?id=3>